

10/774, 082  
Search  
LYCOK 7/7/07

d his

(FILE 'HOME' ENTERED AT 13:45:08 ON 07 JUL 2007)

FILE 'BIOSIS, CAPLUS, EMBASE, MEDLINE, JAPIO' ENTERED AT 13:45:22 ON 07  
JUL 2007

L1 26 S GRADIFLOW AND ANTIBOD?  
L2 14 DUPLICATE REMOVE L1 (12 DUPLICATES REMOVED)  
L3 0 S L2 AND PD<1998  
L4 140 S GRADIFLOW?  
L5 23 S L4 AND PD<1999  
L6 7 DUPLICATE REMOVE L5 (16 DUPLICATES REMOVED)  
L7 7 S L6 AND PROTEIN?

=>

d his

(FILE 'HOME' ENTERED AT 13:45:08 ON 07 JUL 2007)

FILE 'BIOSIS, CAPLUS, EMBASE, MEDLINE, JAPPIO' ENTERED AT 13:45:22 ON 07 JUL 2007

L1	26 S GRADIFLOW AND ANTIBOD?
L2	14 DUPLICATE REMOVE L1 (12 DUPLICATES REMOVED)
L3	0 S L2 AND PD<1998
L4	140 S GRADIFLOW?
L5	23 S L4 AND PD<1999
L6	7 DUPLICATE REMOVE L5 (16 DUPLICATES REMOVED)
L7	7 S L6 AND PROTEIN?

=>

AN 1997:385617 BIOSIS

DN PREV199799684820

TI Purification by reflux electrophoresis of whey proteins and of a recombinant protein expressed in Dictyostelium discoideum.

AU Corthals, Garry L.; Collins, Brett M.; Mabbutt, Bridget C.; Williams, Keith L. [Reprint author]; Gooley, Andrew A.

CS Macquarie Univ. Centre Analytical Biotechnol., Sch. Biological Sci., Macquarie Univ., Sydney, NSW 2109, Australia

SO Journal of Chromatography A, (1997) Vol. 773, No. 1-2, pp. 299-309.

CODEN: JOCRAM. ISSN: 0021-9673.

DT Article

LA English

ED Entered STN: 10 Sep 1997

Last Updated on STN: 10 Sep 1997

AB Protein purification that combines the use of molecular mass exclusion membranes with electrophoresis is particularly powerful as it uses properties inherent to both techniques. The use of membranes allows efficient processing and is easily scaled up, while electrophoresis permits high resolution separation under mild conditions. The Gradiflow apparatus combines these two technologies as it uses polyacrylamide membranes to influence electrokinetic separations. The reflux electrophoresis process consists of a series of cycles incorporating a forward phase and a reverse phase. The forward phase involves collection of a target protein that passes through a separation membrane before trailing proteins in the same solution. The forward phase is repeated following clearance of the membrane in the reverse phase by reversing the current. We have devised a strategy to establish optimal reflux separation parameters, where membranes are chosen for a particular operating range and protein transfer is monitored at different pH values. In addition, forward and reverse phase times are determined during this process. Two examples of the reflux method are described. In the first case, we describe the purification strategy for proteins from a complex mixture which contains proteins of higher electrophoretic mobility than the target protein. This is a two-step procedure, where first proteins of higher mobility than the target protein are removed from the solution by a series of reflux cycles, so that the target protein remains as the leading fraction. In the second step the target protein is collected, as it has become the leading fraction of the remaining proteins. In the second example we report the development of a reflux strategy which allowed a rapid one-step preparative purification of a recombinant protein, expressed in Dictyostelium discoideum. These strategies demonstrate that the Gradiflow is amenable to a wide range of applications, as the protein of interest is not necessarily required to be the leading fraction in solution.

CC Genetics - Plant 03504

Genetics - Animal 03506

Biochemistry methods - Proteins, peptides and amino acids 10054

Biophysics - Methods and techniques 10504

Food microbiology - Biosynthesis, bioassay and fermentation 39007

Plant physiology - Chemical constituents 51522

Plant physiology - Apparatus and methods 51524

Invertebrata: comparative, experimental morphology, physiology and pathology - Protozoa 64002

IT Major Concepts

Biochemistry and Molecular Biophysics; Bioprocess Engineering;

Genetics; Methods and Techniques; Physiology

IT Miscellaneous Descriptors

BIOCHEMISTRY AND BIOPHYSICS; EXPRESSION SYSTEM; METHODOLOGY;

PURIFICATION; PURIFICATION METHOD; RECOMBINANT PROTEINS;

REFLUX ELECTROPHORESIS; WHEY PROTEINS

✓ pulled  
LYCOK

ORGN Classifier

Myxophyta 15700

Super Taxa

Fungi; Plantae

Organism Name

Dictyostelium discoideum

Myxophyta

Taxa Notes

Fungi, Microorganisms, Nonvascular Plants, Plants

ORGN Classifier

Sarcodina 35300

Super Taxa

Protozoa; Invertebrata; Animalia

Organism Name

Sarcodina

Taxa Notes

Animals, Invertebrates, Microorganisms, Protozoans

ANSWER 1 OF 7 BIOSIS COPYRIGHT (c) 2007 The Thomson Corporation on STN

AN 1997:385617 BIOSIS

DN PREV199799684820

TI Purification by reflux electrophoresis of whey proteins and of a recombinant protein expressed in *Dictyostelium discoideum*.

AU Corthals, Garry L.; Collins, Brett M.; Mabbutt, Bridget C.; Williams, Keith L. [Reprint author]; Gooley, Andrew A.

CS Macquarie Univ. Centre Analytical Biotechnol., Sch. Biological Sci., Macquarie Univ., Sydney, NSW 2109, Australia

SO Journal of Chromatography A, (1997) Vol. 773, No. 1-2, pp. 299-309.  
CODEN: JOCRAM. ISSN: 0021-9673.

DT Article

LA English

ED Entered STN: 10 Sep 1997  
Last Updated on STN: 10 Sep 1997

AB Protein purification that combines the use of molecular mass exclusion membranes with electrophoresis is particularly powerful as it uses properties inherent to both techniques. The use of membranes allows efficient processing and is easily scaled up, while electrophoresis permits high resolution separation under mild conditions. The Gradiflow apparatus combines these two technologies as it uses polyacrylamide membranes to influence electrokinetic separations. The reflux electrophoresis process consists of a series of cycles incorporating a forward phase and a reverse phase. The forward phase involves collection of a target protein that passes through a separation membrane before trailing proteins in the same solution. The forward phase is repeated following clearance of the membrane in the reverse phase by reversing the current. We have devised a strategy to establish optimal reflux separation parameters, where membranes are chosen for a particular operating range and protein transfer is monitored at different pH values. In addition, forward and reverse phase times are determined during this process. Two examples of the reflux method are described. In the first case, we describe the purification strategy for proteins from a complex mixture which contains proteins of higher electrophoretic mobility than the target protein. This is a two-step procedure, where first proteins of higher mobility than the target protein are removed from the solution by a series of reflux cycles, so that the target protein remains as the leading fraction. In the second step the target protein is collected, as it has become the leading fraction of the remaining proteins. In the second example we report the development of a reflux strategy which allowed a rapid one-step preparative purification of a recombinant protein, expressed in *Dictyostelium discoideum*. These strategies demonstrate that the Gradiflow is amenable to a wide range of applications, as the protein of interest is not necessarily required to be the leading fraction in solution.

CC Genetics - Plant 03504  
Genetics - Animal 03506  
Biochemistry methods - Proteins, peptides and amino acids 10054  
Biophysics - Methods and techniques 10504  
Food microbiology - Biosynthesis, bioassay and fermentation 39007  
Plant physiology - Chemical constituents 51522  
Plant physiology - Apparatus and methods 51524  
Invertebrata: comparative, experimental morphology, physiology and pathology - Protozoa 64002

IT Major Concepts  
Biochemistry and Molecular Biophysics; Bioprocess Engineering; Genetics; Methods and Techniques; Physiology

IT Miscellaneous Descriptors  
BIOCHEMISTRY AND BIOPHYSICS; EXPRESSION SYSTEM; METHODOLOGY; PURIFICATION; PURIFICATION METHOD; RECOMBINANT PROTEINS; REFLUX ELECTROPHORESIS; WHEY PROTEINS

ORGN Classifier  
Myxophyta 15700  
Super Taxa  
Fungi; Plantae  
Organism Name  
Dictyostelium discoideum  
Myxophyta  
Taxa Notes  
Fungi, Microorganisms, Nonvascular Plants, Plants  
ORGN Classifier  
Sarcodina 35300  
Super Taxa  
Protozoa; Invertebrata; Animalia  
Organism Name  
Sarcodina  
Taxa Notes  
Animals, Invertebrates, Microorganisms, Protozoans

AN 1996:181832 BIOSIS  
 DN PREV199698737961  
 TI Preparative affinity membrane electrophoresis.  
 AU Horvath, Z. Stephen [Reprint author]; Gooley, Andrew A.; Wrigley, Colin W.; Margolis, Joel; Williams, Keith L.  
 CS Macquarie Univ. Cent. Analytical Biochem., North Ryde, NSW 2113, Australia  
 SO Electrophoresis, (1996) Vol. 17, No. 1, pp. 224-226.  
 CODEN: ELCTDN. ISSN: 0173-0835.  
 DT Article  
 LA English  
 ED Entered STN: 29 Apr 1996  
 Last Updated on STN: 29 Apr 1996  
 AB Using the patented Gradiflow system in conjunction with newly developed affinity membranes, the suitability of an electrokinetic technique for affinity fractionation was investigated. Blue dextran incorporated into an affinity membrane was used to deplete a solution of horse serum of albumin, with the result that a majority of serum proteins were enriched tenfold relative to albumin. The technique, when fully developed, would offer some advantages over affinity chromatography, since to a degree it is possible to control which components of the sample are presented to the affinity matrix. Furthermore, the technique would extend the capabilities of the already multifunctional Gradiflow system.  
 CC Biochemistry methods - Proteins, peptides and amino acids 10054  
 Biochemistry studies - Proteins, peptides and amino acids 10064  
 Biophysics - Methods and techniques 10504  
 Biophysics - Membrane phenomena 10508  
 IT Major Concepts  
 Biochemistry and Molecular Biophysics; Membranes (Cell Biology);  
 Methods and Techniques  
 IT Miscellaneous Descriptors  
 ANALYTICAL METHOD; PROTEIN; PURIFICATION METH

✓ pulled

AN 1996:181832 BIOSIS  
 DN PREV199698737961  
 TI Preparative affinity membrane electrophoresis.  
 AU Horvath, Z. Stephen [Reprint author]; Gooley, Andrew A.; Wrigley, Colin W.; Margolis, Joel; Williams, Keith L.  
 CS Macquarie Univ. Cent. Analytical Biochem., North Ryde, NSW 2113, Australia  
 SO Electrophoresis, (1996) Vol. 17, No. 1, pp. 224-226.  
 CODEN: ELCTDN. ISSN: 0173-0835.  
 DT Article  
 LA English  
 ED Entered STN: 29 Apr 1996  
 Last Updated on STN: 29 Apr 1996  
 AB Using the patented Gradiflow system in conjunction with newly developed affinity membranes, the suitability of an electrokinetic technique for affinity fractionation was investigated. Blue dextran incorporated into an affinity membrane was used to deplete a solution of horse serum of albumin, with the result that a majority of serum proteins were enriched tenfold relative to albumin. The technique, when fully developed, would offer some advantages over affinity chromatography, since to a degree it is possible to control which components of the sample are presented to the affinity matrix. Furthermore, the technique would extend the capabilities of the already multifunctional Gradiflow system.  
 CC Biochemistry methods - Proteins, peptides and amino acids 10054  
 Biochemistry studies - Proteins, peptides and amino acids 10064  
 Biophysics - Methods and techniques 10504  
 Biophysics - Membrane phenomena 10508  
 IT Major Concepts  
 Biochemistry and Molecular Biophysics; Membranes (Cell Biology);  
 Methods and Techniques  
 IT Miscellaneous Descriptors  
 ANALYTICAL METHOD; PROTEIN; PURIFICATION METH



✓ pulled  
WC

ANSWER 5 OF 7 BIOSIS COPYRIGHT (c) 2007 The Thomson Corporation on STN

AN 1995:16043 BIOSIS

DN PREV199598030343

TI Multifunctional apparatus for electrokinetic processing of proteins.

AU Horvath, Z. Stephen; Corthals, Garry L.; Wrigley, Colin W.; Margolis, Joel

CS Macquarie Univ. Cent. Analytical Biotechnol., Sydney NSW 2109, Australia

SO Electrophoresis, (1994) Vol. 15, No. 7, pp. 968-971.

CODEN: ELCTDN. ISSN: 0173-0835.

DT Article

LA English

ED Entered STN: 11 Jan 1995

Last Updated on STN: 11 Jan 1995

CC Biochemistry methods - Nucleic acids, purines and pyrimidines 10052

Biochemistry studies - Proteins, peptides and amino acids 10064

Biophysics - Methods and techniques 10504

Biophysics - Molecular properties and macromolecules 10506

Biophysics - Membrane phenomena 10508

External effects - Electric, magnetic and gravitational phenomena 10610

IT Major Concepts

Biochemistry and Molecular Biophysics; Membranes (Cell Biology);

Methods and Techniques

IT Miscellaneous Descriptors

CHARGE; CONCENTRATION; ELECTRODIALYSIS; GRADIFLOW;

PURIFICATION METHOD; SIZE

ANSWER 5 OF 7 BIOSIS COPYRIGHT (c) 2007 The Thomson Corporation on STN  
AN 1995:16043 BIOSIS  
DN PREV199598030343  
TI Multifunctional apparatus for electrokinetic processing of  
proteins.  
AU Horvath, Z. Stephen; Corthals, Garry L.; Wrigley, Colin W.; Margolis, Joel  
CS Macquarie Univ. Cent. Analytical Biotechnol., Sydney NSW 2109, Australia  
SO Electrophoresis, (1994) Vol. 15, No. 7, pp. 968-971.  
CODEN: ELCTDN. ISSN: 0173-0835.  
DT Article  
LA English.  
ED Entered STN: 11 Jan 1995  
Last Updated on STN: 11 Jan 1995  
CC Biochemistry methods - Nucleic acids, purines and pyrimidines 10052  
Biochemistry studies - Proteins, peptides and amino acids 10064  
Biophysics - Methods and techniques 10504  
Biophysics - Molecular properties and macromolecules 10506  
Biophysics - Membrane phenomena 10508  
External effects - Electric, magnetic and gravitational phenomena 10610  
IT Major Concepts  
Biochemistry and Molecular Biophysics; Membranes (Cell Biology);  
Methods and Techniques  
IT Miscellaneous Descriptors  
CHARGE; CONCENTRATION; ELECTRODIALYSIS; GRADIFLOW;  
PURIFICATION METHOD; SIZE

✓pulled

ANSWER 3 OF 7 BIOSIS COPYRIGHT (c) 2007 The Thomson Corporation on STN  
AN 1996:375530 BIOSIS  
DN PREV199699097886  
TI The role of pH and membrane porosity in preparative electrophoresis.  
AU Corthals, Garry L.; Margolis, Joel; Williams, Keith L. [Reprint author];  
Gooley, Andrew A.  
CS Macquarie Univ. Centre Analytical Biotechnology, Sch. Biological Sci.,  
Macquaire Univ., North Ryde N.S.W., 2109, Sydney, Australia  
SO Electrophoresis, (1996) Vol. 17, No. 4, pp. 771-775.  
CODEN: ELCTDN. ISSN: 0173-0835.  
DT Article  
LA English.  
ED Entered STN: 26 Aug 1996  
Last Updated on STN: 26 Aug 1996  
AB The Gradiflow is a preparative electrophoresis apparatus,  
allowing fractionation based on a combination of size and charge of  
proteins in their native (unreduced) form. The preparative  
fractionation of two proteins of similar size and isoelectric  
point is demonstrated using the Gradiflow. A separation  
membrane of appropriate pore size was chosen and then fractionation was  
"fine tuned" by selecting an appropriate buffer pH to accentuate charge  
differences between the proteins of interest. Complete  
separation of mg quantities of bovine serum albumin and ovalbumin was  
achieved within 40 min.\*  
CC Biochemistry methods - Proteins, peptides and amino acids 10054  
Biochemistry studies - General 10060  
Biochemistry studies - Proteins, peptides and amino acids 10064  
Biophysics - Methods and techniques 10504  
Biophysics - Molecular properties and macromolecules 10506  
Biophysics - Membrane phenomena 10508  
IT Major Concepts  
Biochemistry and Molecular Biophysics; Membranes (Cell Biology)  
IT Miscellaneous Descriptors  
ALBUMIN; ANALYTICAL METHOD; GRADIFLOW; OVALBUMIN

ANSWER 3 OF 7 BIOSIS COPYRIGHT (c) 2007 The Thomson Corporation on STN

AN 1996:375530 BIOSIS

DN PREV199699097886

TI The role of pH and membrane porosity in preparative electrophoresis.

AU Corthals, Garry L.; Margolis, Joel; Williams, Keith L. [Reprint author];  
Gooley, Andrew A.

CS Macquarie Univ. Centre Analytical Biotechnology, Sch. Biological Sci.,  
Macquaire Univ., North Ryde N.S.W., 2109, Sydney, Australia

SO Electrophoresis, (1996) Vol. 17, No. 4, pp. 771-775.  
CODEN: ELCTDN. ISSN: 0173-0835.

DT Article

LA English

ED Entered STN: 26 Aug 1996  
Last Updated on STN: 26 Aug 1996

AB The Gradiflow is a preparative electrophoresis apparatus,  
allowing fractionation based on a combination of size and charge of  
proteins in their native (unreduced) form. The preparative  
fractionation of two proteins of similar size and isoelectric  
point is demonstrated using the Gradiflow. A separation  
membrane of appropriate pore size was chosen and then fractionation was  
"fine tuned" by selecting an appropriate buffer pH to accentuate charge  
differences between the proteins of interest. Complete  
separation of mg quantities of bovine serum albumin and ovalbumin was  
achieved within 40 min.\*

CC Biochemistry methods - Proteins, peptides and amino acids 10054  
Biochemistry studies - General 10060  
Biochemistry studies - Proteins, peptides and amino acids 10064  
Biophysics - Methods and techniques 10504  
Biophysics - Molecular properties and macromolecules 10506  
Biophysics - Membrane phenomena 10508

IT Major Concepts  
Biochemistry and Molecular Biophysics; Membranes (Cell Biology)

IT Miscellaneous Descriptors  
ALBUMIN; ANALYTICAL METHOD; GRADIFLOW; OVALBUMIN

✓ pulled  
4/cool

ANSWER 2 OF 7 BIOSIS COPYRIGHT (c) 2007 The Thomson Corporation on STN

AN 1997:290572 BIOSIS

DN PREV199799589775

TI Prefractionation of protein samples prior to two-dimensional electrophoresis.

AU Corthals, Garry L.; Molloy, Mark P.; Herbert, Ben R.; Williams, Keith L. [Reprint author]; Gooley, Andrew A.

CS MUCAB/APAF, Sch. Biol. Sci., Macquarie Univ., Sydney, NSW 2109, Australia

SO Electrophoresis, (1997) Vol. 18, No. 3-4, pp. 317-323.

CODEN: ELCTDN. ISSN: 0173-0835.

DT Article

LA English

ED Entered STN: 9 Jul 1997

Last Updated on STN: 9 Jul 1997

AB Thousands of proteins may be visualized on a two-dimensional (2-D) gel, but only hundreds are present at levels sufficient for chemical analysis. Therefore, prefractionation of protein samples prior to 2-D polyacrylamide gel electrophoresis (PAGE) will be important for the investigation of proteins that are present at sub-picogram levels in physiological samples. We describe an approach to prefractionate protein samples prior to 2-D PAGE using the Gradiflow, which is a new (preparative) electrokinetic membrane apparatus designed to fractionate proteins in a number of different ways. We have fractionated human serum under nonreducing conditions using the 'reflux' mode, in which proteins are fractionated according to their relative mobility under controlled electrophoretic conditions, where the current is periodically reversed. We describe how fractionation occurs and present examples of enrichment of specific proteins.

CC Biochemistry methods - General 10050

Biochemistry studies - General 10060

IT Major Concepts

Biochemistry and Molecular Biophysics

IT Miscellaneous Descriptors

BIOCHEMISTRY AND BIOPHYSICS; CYTOCHEMICAL METHOD; METHODOLOGY;  
PROTEIN PREFRACTIONATION; PROTEIN SAMPLES; TWO  
DIMENSIONAL ELECTROPHORESIS

ANSWER 2 OF 7 BIOSIS COPYRIGHT (c) 2007 The Thomson Corporation on STN  
AN 1997:290572 BIOSIS  
DN PREV199799589775  
TI Prefractionation of protein samples prior to two-dimensional  
electrophoresis.  
AU Corthals, Garry L.; Molloy, Mark P.; Herbert, Ben R.; Williams, Keith L.  
[Reprint author]; Gooley, Andrew A.  
CS MUCAB/APAF, Sch. Biol. Sci., Macquarie Univ., Sydney, NSW 2109, Australia  
SO Electrophoresis, (1997) Vol. 18, No. 3-4, pp. 317-323.  
CODEN: ELCTDN. ISSN: 0173-0835.  
DT Article  
LA English  
ED Entered STN: 9 Jul 1997  
Last Updated on STN: 9 Jul 1997  
AB Thousands of proteins may be visualized on a two-dimensional  
(2-D) gel, but only hundreds are present at levels sufficient for chemical  
analysis. Therefore, prefractionation of protein samples prior  
to 2-D polyacrylamide gel electrophoresis (PAGE) will be important for the  
investigation of proteins that are present at sub-picogram  
levels in physiological samples. We describe an approach to  
prefractionate protein samples prior to 2-D PAGE using the  
Gradiflow, which is a new (preparative) electrokinetic membrane  
apparatus designed to fractionate proteins in a number of  
different ways. We have fractionated human serum under nonreducing  
conditions using the 'reflux' mode, in which proteins are  
fractionated according to their relative mobility under controlled  
electrophoretic conditions, where the current is periodically reversed.  
We describe how fractionation occurs and present examples of enrichment of  
specific proteins.  
CC Biochemistry methods - General 10050  
Biochemistry studies - General 10060  
IT Major Concepts  
Biochemistry and Molecular Biophysics  
IT Miscellaneous Descriptors  
BIOCHEMISTRY AND BIOPHYSICS; CYTOCHEMICAL METHOD; METHODOLOGY;  
PROTEIN PREFRACTIONATION; PROTEIN SAMPLES; TWO  
DIMENSIONAL ELECTROPHORESIS

AN 1992:343363 BIOSIS  
 DN PREV199294035588; BA94:35588  
 TI RAPID TEN-MINUTE PORE-GRADIENT ELECTROPHORESIS OF PROTEINS AND  
 PEPTIDES IN MICROGRAD GELS.  
 AU WRIGLEY C W [Reprint author]; MARGOLIS J  
 CS CSIRO WHEAT RES UNIT, DIV PLANT INDUSTRY, PO BOX 7, NORTH RYDE, NSW 2113,  
 AUST  
 SO Applied and Theoretical Electrophoresis, (1992) Vol. 3, No. 1,  
 pp. 13-16.  
 CODEN: ATELEM. ISSN: 0954-6642.  
 DT Article  
 FS BA  
 LA ENGLISH  
 ED Entered STN: 29 Jul 1992  
 Last Updated on STN: 29 Jul 1992  
 AB Precast gradient gels of short migration length (25 mm) have been  
 developed to provide rapid electrophoretic separation without loss of  
 resolution. These Micrograd gels have been prepared in gel ranges  
 (conventional and unique) to match pore-gradient electrophoresis  
 conditions to proteins/peptides ranging in size from several  
 hundreds to millions. The Hylinx Micrograd gel combines an extreme gel  
 range (6 to 48% polyacrylamide) with a novel crosslinker to provide  
 sieving of polypeptides, and pore-limit electrophoresis of the smallest  
 proteins (e.g. insulin monomer). All gel ranges (such as 3 to  
 30%) provide zone sharpening in routine analysis of conventional  
 protein mixtures (e.g. serum) within 10 min electrophoresis at 200  
 to 300 volts. The gels are thin (1 mm) and thus stain quickly, but the  
 gel cassette is of conventional overall width (83 mm), thus fitting many  
 apparatus designs and accommodating 12 samples. The gels are finding  
 valuable use in screening applications, requiring the electrophoretic  
 analysis of many samples, and in cases where a rapid answer is needed,  
 such as monitoring protein purification. The gels have proved  
 particularly useful, in-house, for the latter application in developing  
 Gradipore's new large-scale preparative electrophoresis system, the  
 Gradiflow.  
 CC Biochemistry methods - Proteins, peptides and amino acids 10054  
 Biochemistry studies - Proteins, peptides and amino acids 10064  
 Biophysics - Methods and techniques 10504  
 Blood - Blood and lymph studies 15002  
 IT Major Concepts  
 Biochemistry and Molecular Biophysics; Blood and Lymphatics (Transport  
 and Circulation)  
 IT Miscellaneous Descriptors  
 PROTEIN FRACTIONATION SERUM PROTEIN PURIFICATION  
 METHOD

AN 1992:343363 BIOSIS  
 DN PREV199294035588; BA94:35588  
 TI RAPID TEN-MINUTE PORE-GRADIENT ELECTROPHORESIS OF PROTEINS AND  
 PEPTIDES IN MICROGRAD GELS.  
 AU WRIGLEY C W [Reprint author]; MARGOLIS J  
 CS CSIRO WHEAT RES UNIT, DIV PLANT INDUSTRY, PO BOX 7, NORTH RYDE, NSW 2113,  
 AUST  
 SO Applied and Theoretical Electrophoresis, (1992) Vol. 3, No. 1,  
 pp. 13-16.  
 CODEN: ATELEM. ISSN: 0954-6642.  
 DT Article  
 FS BA  
 LA ENGLISH  
 ED Entered STN: 29 Jul 1992  
 Last Updated on STN: 29 Jul 1992  
 AB Precast gradient gels of short migration length (25 mm) have been  
 developed to provide rapid electrophoretic separation without loss of  
 resolution. These Micrograd gels have been prepared in gel ranges  
 (conventional and unique) to match pore-gradient electrophoresis  
 conditions to proteins/peptides ranging in size from several  
 hundreds to millions. The Hylinx Micrograd gel combines an extreme gel  
 range (6 to 48% polyacrylamide) with a novel crosslinker to provide  
 sieving of polypeptides, and pore-limit electrophoresis of the smallest  
 proteins (e.g. insulin monomer). All gel ranges (such as 3 to  
 30%) provide zone sharpening in routine analysis of conventional  
 protein mixtures (e.g. serum) within 10 min electrophoresis at 200  
 to 300 volts. The gels are thin (1 mm) and thus stain quickly, but the  
 gel cassette is of conventional overall width (83 mm), thus fitting many  
 apparatus designs and accommodating 12 samples. The gels are finding  
 valuable use in screening applications, requiring the electrophoretic  
 analysis of many samples, and in cases where a rapid answer is needed,  
 such as monitoring protein purification. The gels have proved  
 particularly useful, in-house, for the latter application in developing  
 Gradipore's new large-scale preparative electrophoresis system, the  
 Gradiflow.  
 CC Biochemistry methods - Proteins, peptides and amino acids 10054  
 Biochemistry studies - Proteins, peptides and amino acids 10064  
 Biophysics - Methods and techniques 10504  
 Blood - Blood and lymph studies 15002  
 IT Major Concepts  
 Biochemistry and Molecular Biophysics; Blood and Lymphatics (Transport  
 and Circulation)  
 IT Miscellaneous Descriptors  
 PROTEIN FRACTIONATION SERUM PROTEIN PURIFICATION  
 METHOD



ANSWER 7 OF 7 CAPLUS COPYRIGHT 2007 ACS on STN

AN 1996:435344 CAPLUS

DN 125:136744

ED Entered STN: 24 Jul 1996

TI Large-scale preparative electrophoresis of proteins

AU Wrigley, Colin W.; Margolis, Joel; Manus, H. Perry

CS Div. Plant Industry, CSIRO, North Ryde, 2113, Australia

SO American Biotechnology Laboratory (1996), 14(6), 8, 12

CODEN: ABLAEY; ISSN: 0749-3223

PB International Scientific Communications

DT Journal

LA English

CC 9-1 (Biochemical Methods)

Section cross-reference(s): 48, 66

AB The LM-1000 Gradiflow system (Gradipore Ltd., Sydney, Australia)

is described for the title process that has the following advantages: it permits purification of 1 species of macromol. from a complex mixture in quantity; it provides high purity for the target macromol.; it gives rapid throughput in the sample ranges of many milligrams to many grams with potential for even larger throughput, preferably on a continuous basis; it works under conditions that would not cause protein denaturation or contamination with extraneous compds.; and it is adaptable to other applications such as salt removal and concentrating functions.

ST protein preparative membrane electrophoresis LM1000  
Gradiflow

IT Concentrators

Membrane, biological

(large-scale preparative electrophoresis of proteins)

IT Proteins, preparation

RL: PUR (Purification or recovery); PREP (Preparation)

(large-scale preparative electrophoresis of proteins)

IT Salts, processes

RL: REM (Removal or disposal); PROC (Process)

(large-scale preparative electrophoresis of proteins)

IT Electrophoresis and Ionophoresis

(preparative, apparatus, large-scale preparative electrophoresis of proteins)

ANSWER 7 OF 7 CAPLUS COPYRIGHT 2007 ACS on STN

AN 1996:435344 CAPLUS

DN 125:136744

ED Entered STN: 24 Jul 1996

TI Large-scale preparative electrophoresis of proteins

AU Wrigley, Colin W.; Margolis, Joel; Manus, H. Perry

CS Div. Plant Industry, CSIRO, North Ryde, 2113, Australia

SO American Biotechnology Laboratory (1996), 14(6), 8, 12

CODEN: ABLAEY; ISSN: 0749-3223

PB International Scientific Communications

DT Journal

LA English

CC 9-1 (Biochemical Methods)

Section cross-reference(s): 48, 66

AB The LM-1000 Gradiflow system (Gradipore Ltd., Sydney, Australia)

is described for the title process that has the following advantages: it permits purification of 1 species of macromol. from a complex mixture in quantity; it provides high purity for the target macromol.; it gives rapid throughput in the sample ranges of many milligrams to many grams with potential for even larger throughput, preferably on a continuous basis; it works under conditions that would not cause protein denaturation or contamination with extraneous compds.; and it is adaptable to other applications such as salt removal and concentrating functions.

ST protein preparative membrane electrophoresis LM1000  
Gradiflow

IT Concentrators

Membrane, biological

(large-scale preparative electrophoresis of proteins)

IT Proteins, preparation

RL: PUR (Purification or recovery); PREP (Preparation)

(large-scale preparative electrophoresis of proteins)

IT Salts, processes

RL: REM (Removal or disposal); PROC (Process)

(large-scale preparative electrophoresis of proteins)

IT Electrophoresis and Ionophoresis

(preparative, apparatus, large-scale preparative electrophoresis of proteins)

10/774,082  
LYCOOLK

d his

(FILE 'HOME' ENTERED AT 16:16:50 ON 07 JUL 2007)

FILE 'BIOSIS, CAPLUS, EMBASE, MEDLINE, JAPIO' ENTERED AT 16:17:14 ON 07 JUL 2007

L1 726 S (ANTIBOD? SEPARAT?)  
L2 13175 S (PROTEIN SEPARAT?)  
L3 17 S L1 AND L2  
L4 4 S L3 AND ELECTRO?  
L5 4 DUPLICATE REMOVE L4 (0 DUPLICATES REMOVED)  
L6 26 S GRADIFLOW AND ANTIBOD?  
L7 14 DUPLICATE REMOVE L6 (12 DUPLICATES REMOVED)  
L8 0 S L7 AND PD<1998  
L9 14 DUPLICATE REMOVE L7 (0 DUPLICATES REMOVED)  
L10 13 S MONOCLONAL AND ASCITIS  
L11 11 DUPLICATE REMOVE L10 (2 DUPLICATES REMOVED)  
L12 7 S L11 AND PD<1998  
L13 79 S L1 AND ELECTROPHOR?  
L14 59 DUPLICATE REMOVE L13 (20 DUPLICATES REMOVED)  
L15 45 S L14 AND PD<1998  
L16 21 S L15 AND PROTEIN?

=>

d his

(FILE 'HOME' ENTERED AT 16:16:50 ON 07 JUL 2007)

FILE 'BIOSIS, CAPLUS, EMBASE, MEDLINE, JAPIO' ENTERED AT 16:17:14 ON 07 JUL 2007

L1	726 S (ANTIBOD? SEPARAT?)
L2	13175 S (PROTEIN SEPARAT?)
L3	17 S L1 AND L2
L4	4 S L3 AND ELECTRO?
L5	4 DUPLICATE REMOVE L4 (0 DUPLICATES REMOVED)
L6	26 S GRADIFLOW AND ANTIBOD?
L7	14 DUPLICATE REMOVE L6 (12 DUPLICATES REMOVED)
L8	0 S L7 AND PD<1998
L9	14 DUPLICATE REMOVE L7 (0 DUPLICATES REMOVED)
L10	13 S MONOCLONAL AND ASCITIS
L11	11 DUPLICATE REMOVE L10 (2 DUPLICATES REMOVED)
L12	7 S L11 AND PD<1998
L13	79 S L1 AND ELECTROPHOR?
L14	59 DUPLICATE REMOVE L13 (20 DUPLICATES REMOVED)
L15	45 S L14 AND PD<1998
L16	21 S L15 AND PROTEIN?

=>

AN 1998:291202 CAPLUS  
 DN 129:119844  
 ED Entered STN: 20 May 1998  
 TI Separation of monoclonal antibody by protein A hollow fiber  
 affinity membrane adsorption  
 AU Bao, Shixiang; Su, Zhiguo  
 CS Department of Bioengineering, South China University of Technology,  
 Canton, 510641, Peop. Rep. China  
 SO Mo Kexue Yu Jishu (1997), 17(4), 21-24  
 CODEN: MKYJEF; ISSN: 0254-6140  
 PB Mo Kexue Yu Jishu Bianjibu  
 DT Journal  
 LA Chinese  
 CC 9-16 (Biochemical Methods)  
 AB The adsorption properties of protein A hollow fiber affinity  
 membrane were studied using human  $\gamma$ -Ig as a model protein.  
 Anti-human chorionic gonadotropin monoclonal antibody was purified from  
 mouse ascites. SDS-polyacrylamide gel electrophoretic anal.  
 showed that the purified antibody has high purity.  
 ST affinity membrane chorionic gonadotropin antibody; hollow fiber affinity  
 membrane antibody sepn  
 IT Proteins, specific or class  
 RL: BSU (Biological study, unclassified); BUU (Biological use,  
 unclassified); BIOL (Biological study); USES (Uses)  
 (A; separation of monoclonal antibody by protein A hollow fiber  
 affinity membrane adsorption)  
 IT Antibodies  
 RL: PEP (Physical, engineering or chemical process); THU (Therapeutic  
 use); BIOL (Biological study); PROC (Process); USES (Uses)  
 (monoclonal; separation of monoclonal antibody by protein A hollow  
 fiber affinity membrane adsorption)  
 IT Membranes, nonbiological  
 Pregnancy  
 (separation of monoclonal antibody by protein A hollow fiber  
 affinity membrane adsorption)  
 IT 9002-61-3, Chorionic gonadotropin  
 RL: BSU (Biological study, unclassified); BIOL (Biological study)  
 (separation of monoclonal antibody by protein A hollow fiber  
 affinity membrane adsorption)

AN 1998:291202 CAPLUS  
 DN 129:119844  
 ED Entered STN: 20 May 1998  
 TI Separation of monoclonal antibody by protein A hollow fiber  
 affinity membrane adsorption  
 AU Bao, Shixiang; Su, Zhiguo  
 CS Department of Bioengineering, South China University of Technology,  
 Canton, 510641, Peop. Rep. China  
 SO Mo Kexue Yu Jishu (1997), 17(4), 21-24  
 CODEN: MKYJEF; ISSN: 0254-6140  
 PB Mo Kexue Yu Jishu Bianjibu  
 DT Journal  
 LA Chinese  
 CC 9-16 (Biochemical Methods)  
 AB The adsorption properties of protein A hollow fiber affinity  
 membrane were studied using human  $\gamma$ -Ig as a model protein.  
 Anti-human chorionic gonadotropin monoclonal antibody was purified from  
 mouse ascites. SDS-polyacrylamide gel electrophoretic anal.  
 showed that the purified antibody has high purity.  
 ST affinity membrane chorionic gonadotropin antibody; hollow fiber affinity  
 membrane antibody sepn  
 IT Proteins, specific or class  
 RL: BSU (Biological study, unclassified); BUU (Biological use,  
 unclassified); BIOL (Biological study); USES (Uses)  
 (A; separation of monoclonal antibody by protein A hollow fiber  
 affinity membrane adsorption)  
 IT Antibodies  
 RL: PEP (Physical, engineering or chemical process); THU (Therapeutic  
 use); BIOL (Biological study); PROC (Process); USES (Uses)  
 (monoclonal; separation of monoclonal antibody by protein A hollow  
 fiber affinity membrane adsorption)  
 IT Membranes, nonbiological  
 Pregnancy  
 (separation of monoclonal antibody by protein A hollow fiber  
 affinity membrane adsorption)  
 IT 9002-61-3, Chorionic gonadotropin  
 RL: BSU (Biological study, unclassified); BIOL (Biological study)  
 (separation of monoclonal antibody by protein A hollow fiber  
 affinity membrane adsorption)

ANSWER 16 OF 21 CAPLUS COPYRIGHT 2007 ACS on STN

AN 1966:450690 CAPLUS

DN 65:50690

OREF 65:9509e-f

ED Entered STN: 22 Apr 2001

TI Preparative electrophoretic separation of rabbit serum proteins and antibodies

AU Freeman, M. J.; Stavitsky, A. B.

CS Western Reserve Univ., Cleveland, OH

SO Immunochemistry (1966), 3(4), 257-66

CODEN: IMCHAZ; ISSN: 0019-2791

DT Journal

LA English

CC 67 (Immunochemistry)

AB The application of a medium-free, continuous-flow electrophoretic separation was evaluated for the characterization and purification of rabbit serum proteins and antibody. Electrophoretic sepns. of large vols. of whole antiserum, serum, or fractions which corresponded well to the distribution of serum proteins in agar gel or other types of electrophoresis were obtained. Considerable electrophoretic heterogeneity was observed for the antibody of 2 pooled sera. The electrophoretic method described is a useful addnl. method for the preliminary physicochem. purification of antibody. 23 references.

IT Proteins

(blood serum, separation by continuous-flow electrophoresis)

IT Antibodies

(sepn. of, by continuous-flow electrophoresis)

ANSWER 16 OF 21 CAPLUS COPYRIGHT 2007 ACS on STN

AN 1966:450690 CAPLUS

DN 65:50690

OREF 65:9509e-f

ED Entered STN: 22 Apr 2001

TI Preparative electrophoretic separation of rabbit serum proteins and antibodies

AU Freeman, M. J.; Stavitsky, A. B.

CS Western Reserve Univ., Cleveland, OH

SO Immunochemistry (1966), 3(4), 257-66

CODEN: IMCHAZ; ISSN: 0019-2791

DT Journal

LA English

CC 67 (Immunochemistry)

AB The application of a medium-free, continuous-flow electrophoretic separation was evaluated for the characterization and purification of rabbit serum proteins and antibody. Electrophoretic sepn. of large vols. of whole antiserum, serum, or fractions which corresponded well to the distribution of serum proteins in agar gel or other types of electrophoresis were obtained. Considerable electrophoretic heterogeneity was observed for the antibody of 2 pooled sera. The electrophoretic method described is a useful addnl. method for the preliminary physicochem. purification of antibody. 23 references.

IT Proteins

(blood serum, separation by continuous-flow electrophoresis)

IT Antibodies

(sepn. of, by continuous-flow electrophoresis)



AN 1990:135559 CAPLUS  
 DN 112:135559  
 ED Entered STN: 13 Apr 1990  
 TI Method and apparatus for continuous isoelectric separation of proteins  
 IN Stimpson, Donald Irvine  
 PA Monsanto Co., USA  
 SO Eur. Pat. Appl., 15 pp.  
 CODEN: EPXXDW  
 DT Patent  
 LA English  
 IC ICM C07K003-16  
 ICS B01D013-00  
 CC 9-1 (Biochemical Methods)  
 FAN.CNT 1

	PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
PI	EP 323948	A2	19890712	EP 1989-870001	19890104 <--
	EP 323948	A3	19911121		
	R: AT, BE, CH, DE, ES, FR, GB, GR, IT, LI, LU, NL, SE				
	US 5114555	A	19920519	US 1988-273780	19881123 <--
	DK 8900019	A	19890706	DK 1989-19	19890104 <--
	NO 8900038	A	19890706	NO 1989-38	19890104 <--
	JP 02006739	A	19900110	JP 1989-45	19890104 <--
PRAI	US 1988-140855	A	19880105		
	US 1988-273780	A	19881123		

## CLASS

PATENT NO.	CLASS	PATENT FAMILY CLASSIFICATION CODES
EP 323948	ICM	C07K003-16
	ICS	B01D013-00
	IPCI	C07K0003-16 [ICM,4]; B01D0013-00 [ICS,4]
	IPCR	G01N0027-447 [I,C*]; G01N0027-447 [I,A]; B01D0057-02 [I,C*]; B01D0057-02 [I,A]; C07K0001-00 [I,C*]; C07K0001-28 [I,A]; C07K0016-00 [I,C*]; C07K0016-00 [I,A]
US 5114555	IPCI	G01N0027-26 [ICM,5]; B01D0057-02 [ICS,5]
	IPCR	G01N0027-447 [I,C*]; G01N0027-447 [I,A]; B01D0057-02 [I,C*]; B01D0057-02 [I,A]; C07K0001-00 [I,C*]; C07K0001-28 [I,A]; C07K0016-00 [I,C*]; C07K0016-00 [I,A]
DK 8900019	NCL	204/601.000; 204/610.000; 204/644.000
	IPCI	B01D0013-01 [ICM,4]
	IPCR	G01N0027-447 [I,C*]; G01N0027-447 [I,A]; B01D0057-02 [I,C*]; B01D0057-02 [I,A]; C07K0001-00 [I,C*]; C07K0001-28 [I,A]; C07K0016-00 [I,C*]; C07K0016-00 [I,A]
NO 8900038	IPCI	C07K0003-14 [ICM,4]; C07K0003-12 [ICS,4]
	IPCR	G01N0027-447 [I,C*]; G01N0027-447 [I,A]; B01D0057-02 [I,C*]; B01D0057-02 [I,A]; C07K0001-00 [I,C*]; C07K0001-28 [I,A]; C07K0016-00 [I,C*]; C07K0016-00 [I,A]
JP 02006739	IPCI	G01N0027-26 [ICM,4]

AB The title method and apparatus, for the continuous separation of a target protein or protein fraction from a protein mixture containing  $\geq 2$  proteins at a pH equal to the pI of the target protein, is provided. The apparatus includes  $\geq 1$  nonionic, nonelec. conductive porous membrane conduit through which the protein mixture-containing solution is passed. The conduit, of polysulfone, polyether sulfone, polypropylene, or polyvinylidene difluoride, is positioned to serve as a septum between the buffer chambers and is adapted to permit free flow of electrophoretically driven proteins across the diameter of the conduit. The membrane is

## ANSWER 10 OF 21 CAPLUS COPYRIGHT 2007 ACS on STN

AN 1990:135559 CAPLUS  
 DN 112:135559  
 ED Entered STN: 13 Apr 1990  
 TI Method and apparatus for continuous isoelectric separation of  
 proteins  
 IN Stimpson, Donald Irvine  
 PA Monsanto Co., USA  
 SO Eur. Pat. Appl., 15 pp.  
 CODEN: EPXXDW  
 DT Patent  
 LA English  
 IC ICM C07K003-16  
 ICS B01D013-00  
 CC 9-1 (Biochemical Methods)  
 FAN.CNT 1

	PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
PI	EP 323948	A2	19890712	EP 1989-870001	19890104 <--
	EP 323948	A3	19911121		
	R: AT, BE, CH, DE, ES, FR, GB, GR, IT, LI, LU, NL, SE				
	US 5114555	A	19920519	US 1988-273780	19881123 <--
	DK 8900019	A	19890706	DK 1989-19	19890104 <--
	NO 8900038	A	19890706	NO 1989-38	19890104 <--
	JP 02006739	A	19900110	JP 1989-45	19890104 <--
PRAI	US 1988-140855	A	19880105		
	US 1988-273780	A	19881123		

## CLASS

PATENT NO.	CLASS	PATENT FAMILY CLASSIFICATION CODES
EP 323948	ICM	C07K003-16
	ICS	B01D013-00
	IPCI	C07K0003-16 [ICM,4]; B01D0013-00 [ICS,4]
	IPCR	G01N0027-447 [I,C*]; G01N0027-447 [I,A]; B01D0057-02 [I,C*]; B01D0057-02 [I,A]; C07K0001-00 [I,C*]; C07K0001-28 [I,A]; C07K0016-00 [I,C*]; C07K0016-00 [I,A]
US 5114555	IPCI	G01N0027-26 [ICM,5]; B01D0057-02 [ICS,5]
	IPCR	G01N0027-447 [I,C*]; G01N0027-447 [I,A]; B01D0057-02 [I,C*]; B01D0057-02 [I,A]; C07K0001-00 [I,C*]; C07K0001-28 [I,A]; C07K0016-00 [I,C*]; C07K0016-00 [I,A]
DK 8900019	NCL	204/601.000; 204/610.000; 204/644.000
	IPCI	B01D0013-01 [ICM,4]
	IPCR	G01N0027-447 [I,C*]; G01N0027-447 [I,A]; B01D0057-02 [I,C*]; B01D0057-02 [I,A]; C07K0001-00 [I,C*]; C07K0001-28 [I,A]; C07K0016-00 [I,C*]; C07K0016-00 [I,A]
NO 8900038	IPCI	C07K0003-14 [ICM,4]; C07K0003-12 [ICS,4]
	IPCR	G01N0027-447 [I,C*]; G01N0027-447 [I,A]; B01D0057-02 [I,C*]; B01D0057-02 [I,A]; C07K0001-00 [I,C*]; C07K0001-28 [I,A]; C07K0016-00 [I,C*]; C07K0016-00 [I,A]
JP 02006739	IPCI	G01N0027-26 [ICM,4]

AB The title method and apparatus, for the continuous separation of a target protein or protein fraction from a protein mixture containing  $\geq 2$  proteins at a pH equal to the pI of the target protein, is provided. The apparatus includes  $\geq 1$  nonionic, nonelec. conductive porous membrane conduit through which the protein mixture-containing solution is passed. The conduit, of polysulfone, polyether sulfone, polypropylene, or polyvinylidene difluoride, is positioned to serve as a septum between the buffer chambers and is adapted to permit free flow of electrophoretically driven proteins across the diameter of the conduit. The membrane is

impregnated with a hydrogel to decrease its hydraulic permeability. The conduit is subjected to the influence of an elec. field substantially perpendicular to fluid flow through it, resulting in movement of all charged proteins from the conduit lumen. The target protein, which is unaffected by the elec. field when the pH is its pI, is collected from the conduit outlet in substantially purified form. Thus, a polypropylene fiber membrane (0.2  $\mu$ m pore diameter) was impregnated with 1% (weight/volume) agarose and the dry membrane was wetted 1st in Me<sub>2</sub>CHOH, then in water. The rewetted membrane was submerged in 1% agarose at .apprx.60° for 16-24 h, then drained of excess hydrogel, cooled, and mounted in the apparatus chamber (schematic diagrams included). Using 0.06 M barbital buffer (pH 8.6) as carrier buffer, the apparatus was used to sep.  $\gamma$ -globulins from a protein mixture also containing cytochrome C and phycocyanin. The presence of  $\gamma$ -globulins inside the membrane and the selective removal of cytochrome C and phycocyanin into the surrounding buffer chambers was confirmed by gel electrophoresis anal.

- ST app continuous isoelec protein sepn; membrane conduit continuous isoelec protein sepn; polypropylene agarose membrane isoelec protein sepn
- IT Proteins, analysis  
RL: ANST (Analytical study)  
(apparatus for continuous isoelec. separation of, pH in relation to)
- IT Membranes  
(as conduit in apparatus for continuous isoelec. protein separation)
- IT Ceramic materials and wares  
Polysulfones, uses and miscellaneous  
RL: USES (Uses)  
(membrane conduit of, in apparatus for continuous isoelec)

impregnated with a hydrogel to decrease its hydraulic permeability. The conduit is subjected to the influence of an elec. field substantially perpendicular to fluid flow through it, resulting in movement of all charged proteins from the conduit lumen. The target protein, which is unaffected by the elec. field when the pH is its pI, is collected from the conduit outlet in substantially purified form. Thus, a polypropylene fiber membrane (0.2  $\mu$ m pore diameter) was impregnated with 1% (weight/volume) agarose and the dry membrane was wetted 1st in Me<sub>2</sub>CHOH, then in water. The rewetted membrane was submerged in 1% agarose at .apprx.60° for 16-24 h, then drained of excess hydrogel, cooled, and mounted in the apparatus chamber (schematic diagrams included). Using 0.06 M barbital buffer (pH 8.6) as carrier buffer, the apparatus was used to sep.  $\gamma$ -globulins from a protein mixture also containing cytochrome C and phycocyanin. The presence of  $\gamma$ -globulins inside the membrane and the selective removal of cytochrome C and phycocyanin into the surrounding buffer chambers was confirmed by gel electrophoresis anal.

- ST app continuous isoelec protein sepn; membrane conduit continuous isoelec protein sepn; polypropylene agarose membrane isoelec protein sepn
- IT Proteins, analysis  
RL: ANST (Analytical study)  
(apparatus for continuous isoelec. separation of, pH in relation to)
- IT Membranes  
(as conduit in apparatus for continuous isoelec. protein separation)
- IT Ceramic materials and wares  
Polysulfones, uses and miscellaneous  
RL: USES (Uses)  
(membrane conduit of, in apparatus for continuous isoelec)